



MARKETING ASSETS

RESTORATIVE KIT BY DR. VON SONTAGH

MEDIA FOLDER
LINK

RESTORATIVE KIT BY DR. VON SONTAGH

MARKETING ASSETS

The Dr. Von Sontagh Restoration Kit contains 4 instruments to perform all aesthetic restorations with precision.

It is available with the CLEANext sterilisation cassette for 5 instruments.

OP80C+ADEP VC, CPRO5C+ADEP O,
CK10+ADEP J, OP30C+ADEP RS



OP80C+ADEP VC | Customized plugger with depth marking

Plugger with depth markings 3 and 5 on the smaller side. \varnothing 1.1/1.8 mm



Round plugger \varnothing 1.8



**\varnothing 1.1
3 and 5 markings**



CPR05C+ADEP O | Customized sharp spatula with ruler

With graduation 1-30 on one side of the 8-panel handle, this graduation helps you to realise a precisely symmetrical smile



1.5mm



1-30 graduation




1.5mm





CK10+ADEP J | Modeling instrument

Consisting of an explorer like tip to model the fissures and place the adhesive with a “like one-hair brush” and a rounded cone on the other side for modelling the slopes.



OP30C+ADEP RS | Multi-purpose sharp spatula

Multi-purpose sharp spatula. It allows spreading, trimming, slicing the composite material. The spatula is very thin with 0.2mm thickness at the tip and slightly curved to follow the anatomy of the tooth and easily shape interdental areas.



RESTORATIVE KIT BY DR. VON SONTAGH

MARKETING ASSETS

 DEPPELER

VIDEO



WATCH THE VIDEO HERE
[HTTPS://YOUTU.BE/EX1VHGD4T3I](https://youtu.be/ex1vhgd4t3i)